

Advances in confocal laser endomicroscopy for neuro-oncological tumor resection: a mini-review of current literature

By Carl Zeiss Meditec AG, Medical Technology

Please note: this brief overview does not substitute reading of the original studies cited, but provides an abridged insight that does not necessarily reflect the opinion of notified bodies and/or regulatory authorities. Please consult the original literature in the references list for full information.

Precise identification of tumor areas is crucial in neuro-oncological tumor resection, particularly in the case of infiltrative gliomas.¹

Conventional neuropathological assessments during neurosurgery are typically based on rapid frozen-section histopathology of small extracted tissue samples using conventional benchtop microscopes. Different staining methods (e.g. hematoxylin and eosin [H&E] staining) or fluorescent dyes are used to provide sufficient image contrast.

A number of neurosurgical studies have indicated a range of benefits of fluorescence-guided microscopy for glial tumor resection, including increased gross total resection rate and progression-free survival.¹² Some studies have even suggested possible increases in overall patient survival.¹² However, fluorescence microscopy does not provide sufficient spatial resolution or sensitivity to effectively visualize the surgical margins of diffuse tumors.

Intraoperative optical-sectioning with handheld probe-based confocal laser endomicroscopy systems provides superior resolution and sensitivity for the detection of infiltrating cells at tumor margins, and can potentially be used to quantify tumor parameters in localized regions of brain tissue during the final stages of tumor resection.¹³

From the viewpoint of neuropathologists, in vivo CLE may offer several benefits as an adjunct to physical histopathological assessment. First, biopsy quality can be optimized through more accurately targeting diseased tissue. Second, imaging can complement physical biopsy investigations by providing information on the architectural context of the tissue, without visual artifacts that are commonly seen in frozen biopsy tissue.⁴ Third, this technology can enable real-time in vivo intraoperative neuropathological consultation, providing practical benefits for tumor resection workflow.^{1,2,12}

Technology

In CLE a low-power laser is used to selectively visualize the area of interest in a tissue on a specific focal plane, providing an 'optical section'.^{3,15} As with a conventional, benchtop confocal microscope, newer systems adapted for surgical endomicroscopy incorporate a pinhole that eliminates light from outside the desired focal plane (background fluorescence), which enables high spatial resolution during 3-dimensional assessments (in X,Y, and Z), as well as the evaluation of tissue architecture at a cellular level.^{2,15}

CLE requires the administration of fluorescent contrast agents immediately before the surgical procedure. Intravenous FNa is generally the preferred fluorophore given its 50-year history of clinical use and a well established safety profile.^{15,16} FNa only penetrates into the CNS in areas where the blood-brain-barrier is compromised, as in most gliomas.^{1,3} Studies have confirmed that pre-operative intravenous FNa administration allows assessment of tissue microstructure (e.g. visualization of tumor cell nuclear and cell dimensions), differentiation of tumor cells from healthy brain tissue, determination of tumor border regions, and visualization of injury to normal brain tissue and brain vasculature.^{1,3}

A number of CLE systems are available, with the confocal tip either integrated into an endoscope (e.g. EC3870K by Pentax, Japan) or into a separate probe (e.g. Cellvizio® by Mauna Kea Technologies, France; ZEISS CONVIVO by Carl Zeiss Meditec AG, Germany). Image data from these systems can be evaluated immediately in the OR, or can be streamed for real-time external evaluation to allow remote neuropathologist support.^{3,4} Thus, unlike current standard surgical practice, CLE allows real-time *in vivo* imaging directly in the relevant tissue, supplementing or potentially even replacing physical biopsies.

Currently CLE technology is mainly used in intra-gastrointestinal diagnosis and polyp removal during routine colonoscopy.¹⁷ However, experience in neurosurgical applications in the brain is increasing.^{1,18}

Animal model studies

Proof-of-concept studies in a range of animal species (mice, rats and pigs) suggest that CLE with real-time neuropathology assessment could increase efficiency in the brain tumor resection workflow, potentially reducing the time required in the OR and associated costs.^{1,4,5} In a laboratory assessment of CLE diagnostic accuracy in a mouse glioma model, neuropathologists were able to differentiate tumor from nontumor tissue with a mean accuracy, specificity, and sensitivity of 90%, 86%, and 96%, respectively, and with high overall interobserver agreement.⁶

Preclinical feasibility studies and initial clinical testing suggest that CLE may enable rapid identification of characteristic (diagnostic) histological features of tumor tissue.^{3,7,18-20}

Human *ex vivo* and *in vivo* studies

A number of human clinical studies including *ex vivo* and *in vivo* evaluations have assessed the utility of CLE for real-time intraoperative differential diagnosis.

The effectiveness of CLE was evaluated in a study based on *ex vivo* human brain biopsies and sections from a number of different tumor types that had distinctive as well as complementary reflectance and fluorescence characteristics. Multimodal imaging including CLE allowed neuropathologists to distinguish gliomas from normal brain tissue and nonglial tumors.²¹ CLE findings were comparable with standard H&E-stained slide assessments, supporting the utility of CLE in differential intraoperative diagnoses.

In a single-center feasibility study with an early version of the current ZEISS CONVIVO CLE system in 31 patients with a variety of brain tumors (mainly meningiomas and metastatic gliomas), the neurosurgeon-neuropathology team were able to identify high-grade gliomas and vascular neo-proliferation and to define tumor margins intraoperatively.⁸ CLE provided sufficient resolution to achieve preliminary diagnosis without removal of biopsy samples.

In a case series based on 50 tumor resections, findings from intraoperative CLE correlated well with corresponding traditional histology in identifying pathognomonic cytoarchitectural features in various brain tumor types; in a blinded sub-analysis, 26 (92.9%) of 28 lesions were diagnosed correctly.⁹

A prospective surgical evaluation of real-time CLE in 74 patients provided histological specificities and sensitivities for gliomas and meningiomas that were comparable to those derived from standard assessments by frozen section histopathology.¹⁰

Supported by findings from preclinical, proof-of-concept data from animal studies,^{5,22} these clinical data suggest that CLE with FNa contrast could allow real-time, interactive, intraoperative identification of brain tumor areas, which may provide improvements in the efficiency of decision making during neuro-oncological tumor resection.^{4,10}

Beyond static 2D imaging

In addition to static two-dimensional image analysis, the Z-stack technology of current CLE systems allows intraoperative three-dimensional construction of time series and video loops to enable assessment of blood cell movements inside the brain vasculature and within brain tissue if oozing blood is present.^{1,2,9} In addition, movement of tumor cells relative to one another and to normal brain cells during intraoperative squeezing can provide a further dimension for neuropathological assessment, by which movement of individual tumor cells in the border regions relative to cells at the tumor core can be assessed.¹

Outlook

Overall, evidence from in vivo and ex vivo clinical studies conducted to date indicates that CLE shows promise for precise identification of tumor margins during assessment and resection procedures.^{1,3} The possibility of real-time in vivo intraoperative neuropathological consultation could bring tangible benefits for the brain tumor surgery workflow,^{1,4,11} and imaging technologies such as CLE could potentially lead to new diagnostic criteria for tumor characterization above and beyond those described using conventional two-dimensional histology.^{5,12} In addition, optical imaging provides information that could potentially be leveraged to reveal previously unrealized features of disease and/or biochemical composition in neurosurgery.²³

Please find further details regarding study design and findings in the **original articles** listed in the **reference list** at the end of this document.

References

The publications below are based on the authors' own professional opinions or their study results. They do not necessarily reflect the opinions of ZEISS and may not be in line with the clinical evaluation or intended purpose of our medical devices. Therefore, suitability of clinical application for each recommendation should be carefully assessed by the concerned physician.

- ¹ Belykh E, Miller EJ, Carotenuto A, Patel AA, Cavallo C, Martirosyan NL, Healey DR, Byvaltsev VA, Scheck AC, Lawton MT, Eschbacher JM, Nakaji P, Preul MC. Progress in confocal laser endomicroscopy for neurosurgery and technical nuances for brain tumor imaging with fluorescein. *Front Oncol* 2019;9:554.
- ² Belykh E, Patel AA, Miller EJ, Bozkurt B, Yağmurlu K, Woolf EC, Scheck AC, Eschbacher JM, Nakaji P, Preul MC. Probe-based three-dimensional confocal laser endomicroscopy of brain tumors: technical note. *Cancer Manag Res* 2018a;10:3109–3123.
- ³ Leierseder S. Confocal endomicroscopy during brain surgery. *Laser+Photonics* 2018. Available at: <https://www.photonik.de/confocal-endomicroscopy-during-brain-surgery/150/21404/363684> [Accessed 2 April 2020].
- ⁴ Fotteler ML, Holl F, Käsbaach C, Schlegel J, Swoboda W. Using confocal endomicroscopy for digital biopsy during brain surgery – presentation of a study protocol. In: *ICT for Health Science Research 2019*, A. Shabo (Shvo) et al. (Eds.), The European Federation for Medical Informatics (EFMI) and IOS Press, Hannover, pp. 237–238.
- ⁵ Sankar T, Delaney PM, Ryan RW, Eschbacher J, Abdelwahab M, Nakaji P, Coons SW, Scheck AC, Smith KA, Spetzler RF, Preul MC. Miniaturized handheld confocal microscopy for neurosurgery: results in an experimental glioblastoma model. *Neurosurgery* 2010;66:410–417.
- ⁶ Belykh E, Miller EJ, Patel AA, Yazdanabadi MI, Martirosyan NL, Yağmurlu K, Bozkurt B, Byvaltsev VA, Eschbacher JM, Nakaji P, Preul MC. Diagnostic accuracy of a confocal laser endomicroscope for in vivo differentiation between normal injured and tumor tissue during fluorescein-guided glioma resection: laboratory investigation. *World Neurosurg* 2018b;115:e337–e348.
- ⁷ Martirosyan NL, Georges J, Eschbacher JM, Cavalcanti DD, Elhadi AM, Abdelwahab MG, Scheck AC, Nakaji P, Spetzler RF, Preul MC. Potential application of a handheld confocal endomicroscope imaging system using a variety of fluorophores in experimental gliomas and normal brain. *Neurosurg Focus* 2014;36:E16.
- ⁸ Sanai N, Eschbacher J, Hattendorf G, Coons SW, Preul MC, Smith KA, Nakaji P, Spetzler RF. Intraoperative confocal microscopy for brain tumors: a feasibility analysis in humans. *Neurosurgery* 2011;68(2 Suppl Operative):282–290.
- ⁹ Eschbacher J, Martirosyan NL, Nakaji P, Sanai N, Preul MC, Smith KA, Coons SW, Spetzler RF. In vivo intraoperative confocal microscopy for real-time histopathological imaging of brain tumors. *J Neurosurg* 2012;116:854–860.
- ¹⁰ Martirosyan NL, Eschbacher JM, Kalani MY, Turner JD, Belykh E, Spetzler RF, Nakaji P, Preul MC. Prospective evaluation of the utility of intraoperative confocal laser endomicroscopy in patients with brain neoplasms using fluorescein sodium: experience with 74 cases. *Neurosurg Focus* 2016;40:E11.
- ¹¹ Belykh E, Cavallo C, Gandhi S, Zhao X, Veljanoski D, Izady Yazdanabadi M, Martirosyan NL, Byvaltsev VA, Eschbacher J, Preul MC, Nakaji P. Utilization of intraoperative confocal laser endomicroscopy in brain tumor surgery. *J Neurosurg Sci* 2018c;62:704–717.
- ¹² Belykh E, Martirosyan NL, Yagmurlu K, Miller EJ, Eschbacher JM, Izady Yazdanabadi M, Bardanova LA, Byvaltsev VA, Nakaji P, Preul MC. Intraoperative fluorescence imaging for personalized brain tumor resection: current state and future directions. *Front Surg* 2016;3:55.

Further references

- ¹³ Wei L, Fujita Y, Sanai N, Liu JTC. Toward quantitative neurosurgical guidance with high-resolution microscopy of 5-aminolevulinic acid-induced protoporphyrin IX. *Front Oncol* 2019;9:592.
- ¹⁴ Hariri LP. In vivo microscopy: will the microscope move from our desk into the patient? *Arch Pathol Lab Med* 2015;139:719–720.
- ¹⁵ Nabi Z, Reddy DN. Optical biopsy in gastroenterology: Focus on confocal laser endomicroscopy. *Indian J Gastroenterol* 2019;38:281–286.
- ¹⁶ Kiesslich R, Burg J, Vieth M, Gnaendiger J, Enders M, Delaney P, Polglase A, McLaren W, Janell D, Thomas S, Nafe B, Galle PR, Neurath MF. Confocal laser endoscopy for diagnosing intraepithelial neoplasias and colorectal cancer in vivo. *Gastroenterology* 2004;127:706–713.
- ¹⁷ Belderbos TDG, van Oijen MGH, Moons LMG, Siersema PD. Implementation of real-time probe-based confocal laser endomicroscopy (pCLE) for differentiation of colorectal polyps during routine colonoscopy. *Endosc Int Open* 2017;5:E1104–E1110.
- ¹⁸ Zehri AH, Ramey W, Georges JF, Mooney MA, Martirosyan NL, Preul MC, Nakaji P. Neurosurgical confocal endomicroscopy: review of contrast agents, confocal systems, and future imaging modalities. *Surg Neurol Int* 2018;5:60.
- ¹⁹ Breuskin D, Divincenzo J, Kim YJ, Urbschat S, Oertel J. Confocal laser endomicroscopy in neurosurgery: a new technique with much potential. *Minim Invasive Surg* 2013;85:1819.
- ²⁰ Krishna SG, Lee JH. Appraisal of needle-based confocal laser endomicroscopy in the diagnosis of pancreatic cysts. *World J Gastroenterol* 2016;22:1701–1710.
- ²¹ Snuderl M, Wirth D, Sheth SA, Bourne SK, Kwon CS, Ancukiewicz M, Curry WT, Frosch MP, Yaroslavsky AN. Dye-enhanced multimodal confocal imaging as a novel approach to intraoperative diagnosis of brain tumors. *Brain Pathol* 2013;23:73–81.
- ²² National Institute of Clinical Excellence (NICE). Cellvizio confocal endomicroscopy system for characterising pancreatic cysts. Medtech innovation briefing [MIB69]. June 2016. Available at: <https://www.nice.org.uk/advice/mib69/chapter/Summary> [Accessed 26 Mar 2020].
- ²³ Wells WA, Thrall M, Sorokina A, Fine J, Krishnamurthy S, Haroon A, Rao B, Shevchuk MM, Wolfsen HC, Tearney GJ, Hariri LP. In vivo and ex vivo microscopy - Moving toward the integration of optical Imaging technologies into pathology practice. *Arch Pathol Lab Med* 2019;143:288–298.

Further references not cited here

Schebesch KM, Rosengarth K, Brawanski A, Proescholdt M, Wendl C, Höhne J, Ott C, Lamecker H, Doenitz C. Clinical benefits of combining different visualization modalities in neurosurgery.

Front Surg 2019;6:56. DOI: 10.3389/fsurg.2019.00056

Izadyyazdanabadi M, Belykh E, Mooney MA, Eschbacher JM, Nakaji P, Yang Y, Preul MC. Prospects for theranostics in neurosurgical imaging: empowering confocal laser endomicroscopy diagnostics via deep learning.

Front Oncol 2018;8:240. DOI: 10.3389/fonc.2018.00240

Peyre M, Clermont-Taranchon E, Stemmer-Rachamimov A, Kalamarides M. Miniaturized handheld confocal microscopy identifies focal brain invasion in a mouse model of aggressive meningioma.

Brain Pathol 2013;23:371–377. DOI: 10.1111/bpa.12039

EN_30_200_00991 CZ:V/2020 International edition: Only for use in selected countries.
Not all products, services or offers are approved or offered in every market and approved labeling and instructions may vary from one country to another.
The statements of the author(s) of this publication reflect only their personal opinion and experience and do not necessarily reflect the opinion of any institution with whom they are affiliated.
The authors alone are responsible for the content of this mini review and any potential resulting infringements. Carl Zeiss Meditec AG does not accept any responsibility or liability of its content.
Published by Carl Zeiss Meditec AG, 2020. www.zeiss.com/med All rights reserved.