

HFA System Software Version 4.1 Offers Important Maintenance Enhancements

The next HFA software, Version 4.1, addresses minor issues noted in the prior Version 4.0 release. One of these issues relates to saving exams on floppy disks. Unlike the previous version software, Version 4.1 enables the user to successfully save exam data to a floppy disk even if a new floppy disk is inserted at the end of a test.

The release notes distributed with the product contain a complete list of software enhancements. All new purchases of the HFA II-*i* perimeter will automatically be shipped with the latest Version 4.1 software already installed in the system. Upgrade kits are now available for HFA II-*i* (700*i* series) users who are currently operating on Software Version 4.0.

Please call the Carl Zeiss Meditec Customer Service Department at 1-800-341-6968 to order your free HFA 4.1 software upgrade kit.



Connect Your HFA II-*i* to Your EMR System with HFA-Net Pro

Building on the network compatible features of HFA-Net, HFA-Net Pro enables bi-directional communication between your HFA II-*i* units and with an increasing number of EMR systems. Connecting to your EMR system will help to increase practice efficiency and enhance workflow by reducing the need for manual data entry and minimizing paper handling.

Most major EMR vendors in the USA are developing interfaces for the HFA — the Gold Standard in perimetry. Ask your EMR vendor if they have an interface available for your HFA, and if not, they may contact Carl Zeiss Meditec for a software development kit. Please note that HFA II models need to be upgraded to the HFA II-*i* series in order to connect to an EMR system with HFA-Net Pro.

The following U.S. companies have completed their HFA interfaces as of December 2006:

Company	Product
OfficeMate Software Solutions (Subsidiary of Marchon Eyewear)	ExamWriter (EMR)
MedFlow	Ophthalmicsuite (EMR) HFA Viewer
Universal Software Solutions (VersaSuite)	VersaSuite (PMS/EMR)

** HFA-Net Pro is available for purchase in the U.S. for HFA II-*i* series perimeters. To check availability of HFA-Net Pro outside of the U.S., please contact your local Carl Zeiss Meditec office.*