

ZEISS AT LISA 809M

ZEISS AT LISA 809MP preloaded

Technical Specifications



Seeing beyond



AT LISA® 809M

Optic Design	Multifocal, diffractive, +3.75 D near add at the IOL plane, aspheric (aberration correcting)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size Diopter Range	1.8 mm
Company Labeled A-Constant ¹	117.8
Diopter Range	0.0 to +32.0 D, 0.5 D increments
ACD	4.75
Refractive Index	1.46
Implantation in	Capsular Bag
Injector/Cartridge Set ²	AT.Shooter A2-2000/AT.Smart Cartridge <i>or</i> AT.Shooter A1-2000/VISCOJECT™-BIO 1.8 Cartridge Set <i>or</i> VISCOJECT-BIO 1.8 Injector Set
Indications	<ul style="list-style-type: none"> ■ Indicated for the visual correction of aphakia secondary to the removal of the crystalline lens in patients with cataract ■ Also indicated for non-cataractous, presbyopic patients who seek greater independence from glasses for intermediate and/or near distances



AT LISA® 809MP preloaded

Optic Design	Multifocal, diffractive, +3.75 D near add at the IOL plane, aspheric (aberration correcting)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size Diopter Range	1.8 mm
Company Labeled A-Constant ¹	117.8
Diopter Range	0.0 to +32.0 D, 0.5 D increments
ACD	4.75
Refractive Index	1.46
Implantation in	Capsular Bag
Injector/Cartridge Set ²	BLUEMIXS 180
Indications	<ul style="list-style-type: none"> ■ Indicated for the visual correction of aphakia secondary to the removal of the crystalline lens in patients with cataract ■ Also indicated for non-cataractous, presbyopic patients who seek greater independence from glasses for intermediate and/or near distances

¹ For optimized A Constants and ACD Constants refer to IOLCon: www.iolcon.org

² Please refer to our web page for the most up-to-date references: www.zeiss.com/injectors



AT LISA 809M/809MP
BLUEMIXS 180 Injector



BLUEMIXS 180 Injector



Carl Zeiss Meditec AG
Goeschwitzer Strasse 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/med/contacts



Carl Zeiss Meditec SAS
27 Avenue Paul Langevin
17180 Périgny
France