



AT TORBI 709M
AT TORBI 709MP preloaded
Technical Specifications



AT TORBI 709M

AT TORBI 709MP preloaded



AT TORBI® 709M

Optic Design	Mono-focal, bitoric, aspheric (aberration neutral)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, bitoric, MICS
Incision Size	1.8 mm up to 24.0 D / 2.2 mm up to 32.0 D
Company Labeled A-Constant ¹	118.3
Diopter Range	Sphere: -10.0 to +32.0 D, 0.5 D increments Cylinder: +1.0 to +12.0 D, 0.5 D increments
ACD	5.14
Implantation in	Bag
Injector / Cartridge Set ²	For IOLs from -10.0 to +28.0 sphere with +1.0 to +4.0 cylinder*: AT.Shooter A2-2000 injector / ACM2 cartridge <i>or</i> AT.Shooter A1-2000 injector / VISCOJECT™-BIO 1.8 Cartridge Set <i>or</i> VISCOJECT-BIO 1.8 Injector Set * For diopters +28.5 to +32.0 D sphere with +1.0 to +12.0 D cylinder use VISCOJECT-BIO 2.2 Injector Set or AT.Shooter A1-2000 injector with VISCOJECT-BIO 2.2 Cartridge Set.

AT TORBI® 709MP preloaded

Optic Design	Mono-focal, bitoric, aspheric (aberration neutral)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface properties
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, bitoric, MICS
Incision Size	1.8 mm
Company Labeled A-Constant ¹	118.3
Diopter Range	Sphere: -10.0 to +28.0 D, 0.5 D increments Cylinder: +1.0 to +4.0 D, 0.5 D increments
ACD	5.14
Implantation in	Bag
Injector / Cartridge Set ²	BLUEMIXS® 180

¹Please refer to our web pages for optimized A-Constants. ²Please refer to our web pages for the most up-to-date references.



0297 AT TORBI 709M / 709MP
ACM2 Cartridge Set



0459 BLUEMIXS 180 injector



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/med/contacts



Carl Zeiss Meditec SAS
27 Avenue Paul Langevin
17180 Périgny - France
Mailing address:
Avenue Paul Langevin, BP5
17053 La Rochelle Cedex 9
France



000000-2009-687



03/2015

EN_32_022_00161

AT TORBI and BLUEMIXS are registered trademarks of Carl Zeiss Meditec. VISCOJECT is a trademark of Medcel (CE 0482). The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper. © 2015 by Carl Zeiss Meditec AG. All copyrights reserved.