

ZEISS AT TORBI 709M/709MP

Monofocal bitoric MICS IOL

AT TORBI® range:

Sph -10.0 to +32.0 D

Cyl +1.0 to +12.0 D

All 0.5 D increments



We make it visible.



ZEISS AT TORBI 709M / 709MP

Fast delivery in 1–2 working days:

Sph +6.0 to +30.0 D




Cyl +1.0 to +6.0 D

Fast delivery for the preloaded AT TORBI 709MP is available Sph +6.0 to +28.0 D and Cyl +1.0 to +4.0 D.

Fast delivery provides next day service to most areas if ordered before the cut-off time. Allow 35 days for toric IOLs outside the fast delivery range. Contact your ZEISS IOL customer service regarding your specific delivery lead times.

// FAST DELIVERY
MADE BY ZEISS

Injector Correlation Chart

Sphere powers (D):	-10.0 to +28.0	+28.5 to +32.0
AT TORBI 709MP preloaded 	BLUEMIXS® 180 injector set Cylinder +1.0 to +4.0 D	Not available
AT TORBI 709M single use injectors 	VISCOJECT™-BIO 1.8 injector and cartridge set Cylinder +1.0 to +4.0 D	VISCOJECT-BIO 2.2 injector and cartridge set* Cylinder +1.0 to +12.0 D
AT TORBI 709M reusable injectors 	AT.Shooter A2-2000 injector with ACM2 cartridge or AT.Shooter A1-2000 injector with VISCOJECT-BIO 1.8 cartridge Cylinder +1.0 to +4.0 D	AT.Shooter A1-2000 injector VISCOJECT-BIO 2.2 cartridge* Cylinder +1.0 to +12.0 D

*Can also be used with spherical power lower than 28.5 D

Z CALC®: precise ZEISS toric IOL calculations.

Go to <https://zcalc.meditec.zeiss.com>

Username:

Password:

Z CALC support email address:

support.zcalc.meditec@zeiss.com



0297

AT TORBI 709M / 709MP
ACM2 Cartridge Set



AT.Shooter A1-2000
AT.Shooter A2-2000



Carl Zeiss Meditec AG

Goeschwitzer Str. 51–52
07745 Jena
Germany

www.meditec.zeiss.com/contacts
www.meditec.zeiss.com/iol



0459

BLUEMIXS 180 injector



Carl Zeiss Meditec SAS

27 Avenue Paul Langevin
17180 Périgny – France

Mailing address:

Avenue Paul Langevin, BP5
17053 La Rochelle Cedex 9
France

www.meditec.zeiss.com/contacts
www.meditec.zeiss.com/iol

EN_32_025_020311

The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. AT TORBI, BLUEMIXS and Z CALC are registered trademarks of Carl Zeiss Meditec. VISCOJECT is a trademark of Medical (CE0482).

© 2014 by Carl Zeiss Meditec AG. All copyrights reserved.



03/2014



000000-2105-383