

# ZEISS CT LUCIA 202

## Technical Specifications



Seeing beyond



### CT LUCIA® 202

Optic Design	Monofocal, spheric
Material	Hydrophobic acrylic
Optic Diameter	6.0 mm
Total Diameter	13.0 mm
Haptic	5°
Lens Design	Three-piece
Incision Size	2.8 mm
Company Labeled A-Constant <sup>1</sup>	118.2
Diopter Range	From +4.0 to +34.0 D, 0.5 D increments
ACD <sup>2</sup>	5.1
Abbe number	50
Implantation in	Capsular bag
Injector/Cartridge Set <sup>2</sup>	DUALTEC™ OD655 Injector for diopter range +4.0 to +34.0 D

<sup>1</sup> For optimized A Constants and ACD Constants refer to IOLCon: [www.iolcon.org](http://www.iolcon.org)

<sup>2</sup> Please refer to our web page for the most up-to-date references: [www.zeiss.com/injectors](http://www.zeiss.com/injectors)

**CE** 0297

CT LUCIA 202



**Carl Zeiss Meditec AG**

Goeschwitzer Strasse 51–52

07745 Jena

Germany

[www.zeiss.com/lucia](http://www.zeiss.com/lucia)

[www.zeiss.com/med/contacts](http://www.zeiss.com/med/contacts)

**EN\_32\_022\_005311** Printed in Germany. CZ-IX/2019 International edition: Only for sale in selected countries.  
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information.  
Subject to changes in design and scope of delivery and due to ongoing technical development. CT LUCIA is either a trademark or registered trademark of Carl Zeiss Meditec AG or other companies of the ZEISS Group in Germany and/or other countries.  
© Carl Zeiss Meditec AG, 2019. All rights reserved.