



CT ASPHINA 604P
Technical Specifications



We make it visible.

CT ASPHINA 604P



CT ASPHINA® 604P

Optic Design	Monofocal, aspheric (aberration correcting)
Material	Hydrophilic acrylic (28 %)
Optic Diameter	6.0 mm
Total Diameter	10.5 mm
Haptic Angulation	0°
Lens Design	Single-piece
Incision Size	2.8 mm
Company Labeled A-Constant ¹	118.1
Diopter Range	From +14.0 to +30.0 D +14.0 to +30.0 D, 0.5 D increments
ACD	5.02
Implantation in	Bag
Injector ²	SkyInvent®

¹Please refer to our web pages for optimized A-Constants. ²Please refer to our web pages for the most up-to-date references.

CE 0459 CT ASPHINA 604P
SkyInvent



Carl Zeiss Meditec AG
Goeschitzer Str. 51–52
07745 Jena
Germany
www.zeiss.com/iol
www.zeiss.com/contacts



Carl Zeiss Meditec SAS
27 Avenue Paul Langevin
17180 Périgny – France
Mailing address:
Avenue Paul Langevin – BP5
17053 La Rochelle – Cedex 9
France

EN_32_022_00951

CT ASPHINA and SkyInvent are registered trademarks of Carl Zeiss Meditec. The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper.
© 2015 by Carl Zeiss Meditec AG. All copyrights reserved.

