



ZEISS PURI CLEAR
Balanced Salt Solution (BSS)
Technical Specifications



ZEISS PURI CLEAR

For intraocular use



PURI CLEAR®	
Application	<p>PURI CLEAR is indicated for the irrigation of the anterior chamber during cataract surgery and other intraocular procedures;</p> <ul style="list-style-type: none">■ for the removal of alpha-chymotrypsin after enzymatic zonulolysis,■ for the irrigation of the anterior chamber after anterior chamber hemorrhage,■ for the irrigation of the cornea during strabotomy and other extraocular procedures,■ for the irrigation of the conjunctiva following application of fluorescein, and■ for moistening the corneal and conjunctival surface during laser treatment.
Composition	<ul style="list-style-type: none">■ Sodium chloride■ Potassium chloride■ Calcium chloride 2-hydrate■ Sodium acetate 3-hydrate■ Trisodium citrate 2-hydrate■ Magnesium chloride 6-hydrate■ in a sterile, isotonic and isosmotic solution <p>The balanced salt solution and its primary packaging are free of latex, phthalates, animal tissues and conservatives.</p>
Physico-chemical data	<p>Osmolality (mOsmol/kg): 295-315 pH value: 6.8-7.6 Refractive index n^o20: 1.3340-1.3360</p>
Properties	<p>Composition, pH, and osmolality of PURI CLEAR correspond to that of the aqueous humor. PURI CLEAR does not contain any active pharmaceutical ingredients.</p>
Package size	<p>Box of:</p> <ul style="list-style-type: none">■ 10 x 250 ml plastic bottle■ 10 x 500 ml plastic bottle■ 10 x 500 ml glass bottle■ 10 x 500 ml infusion bag
Storage temperature (°C)	<p>Plastic bottle: no requirement Glass bottle: 5-25 Infusion bag: no requirement</p>
Shelf life	<p>Plastic bottle: 36 months Glass bottle: 12 months Infusion bag: 36 months</p>

CE 0297

PURI CLEAR



Carl Zeiss Meditec AG

Goeschwitzer Str. 51-52

07745 Jena

Germany

www.zeiss.com/bss

www.zeiss.com/contact

EN_32_022_010311 Printed in Germany. CZ-X/2018 International edition: Only for sale in selected countries.
The contents of the brochure may differ from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information.
Subject to changes in design and scope of delivery and due to ongoing technical development. PURI CLEAR is a registered trademark of Carl Zeiss Meditec AG in Germany and/or other countries.
© Carl Zeiss Meditec AG, 2018. All rights reserved.